

Zeeman Effect Analyzer

Slides made and presented by - Abdallah W. Mahmoud 900221058
- Yousif Shaaban 900232840

Overview

Theoretical Background


- Theory behind the Zeeman effect
- Description and significance of the Bohr Magneton and specific charge

Modules Description

- Description of all the modules of the program
- Significance and functionality of each module


Usage

- Simulation of the code on a sample data
- Verification of the program's results with theoretical results



Theoretical Background

In this section, we explain the Zeeman effect in detail



- The splitting of atomic energy levels, and therefore the spectral lines, of an atom when subjected to an external magnetic field.
- The Zeeman effect provided proof of the existence of the electron's magnetic moment, given by $\mu_J = \mu_s + \mu_L$

$$\mu_L = -\left(\frac{e}{2m_e}\right)L = -\mu_B \left(\frac{L}{\hbar}\right)$$


$$\mu_S = -g_S \left(\frac{e}{2m_e}\right)S = -g_S \mu_B \left(\frac{S}{\hbar}\right)$$

$$\Delta E = -\mu_{\{J,z\}} B \quad \Delta E = g_J \mu_B B M_J$$

$$\mu_{\{J,z\}} = -g_J \mu_B \left(\frac{J_z}{\hbar}\right) = -g_J \mu_B M_J$$

Software Architecture

Flowchart



Main Window

```
1  def __init__(self):
2      super().__init__()
3
4      # Initialize state variables
5      self.images = []
6
6      self.current_image_index = -1
7      self.calibration_distance_mm = 2.0
8      self.measurements = []
9
10     # Initialize components
11     self.ui_manager = UIManager(self)
12     self.ui_manager.setup_layout()
13     self.image_processor = ImageProcessor()
14     self.image_display_manager = ImageDisplayManager(self.image_display,
15             self, self.ui_manager)
15     self.measurement_controller = MeasurementController(self,
16             self.ui_manager)
```

Measurement Control

```
1  def set_mode(self, mode: Optional[str]):  
2      intended_mode = mode  
3  
4      if intended_mode == 'center':  
5          self.initialize_for_new_measurement()  
6          self.current_mode = 'center'  
7      elif intended_mode and intended_mode.startswith('auto_'):  
8          if self.current_measurement.get('center') is None:  
9              QMessageBox.information(self.mw, 'Set Center First',  
10                             'Please set the center point before  
defining an annulus for auto-detection.')  
11         self.current_mode = None  
12         self.is_defining_annulus = False  
13         self.auto_detect_limits = {'lower': None, 'upper': None}  
14         self.mw.update_display()  
15         return  
16     else:  
17         self.current_mode = intended_mode  
18         self.is_defining_annulus = True  
19         self.auto_detect_limits = {'lower': None, 'upper': None}  
20     else:  
21         self.current_mode = intended_mode  
22         self.is_defining_annulus = False  
23         self.auto_detect_limits = {'lower': None, 'upper': None}  
24  
25     self.mw.update_display()
```

Image Processor

Physics Module

```
1  @dataclass
2  class ZeemanMeasurement:
3      B_field: float # Magnetic field strength in Tesla
4
5      wavelength: float # Central wavelength in meters
6      R_center: Optional[float] = None # Radius of central line in mm
7      R_inner: Optional[float] = None # Radius of inner line in mm
8      R_outer: Optional[float] = None # Radius of outer line in mm
9
10     # [Additional fields omitted for brevity]
11
12     def calculate_wavelength_shift(beta: float, beta_c: float, wavelength: float) -> float:
13         """Calculate the wavelength shift based on the refracted angles."""
14         return wavelength * (np.cos(beta_c) / np.cos(beta) - 1)
15
16     def calculate_energy_shift(delta_lambda: float, wavelength: float) -> float:
17         """Calculate the energy shift based on the wavelength shift."""
18         return PLANCK * LIGHT_SPEED * delta_lambda / (wavelength ** 2)
```

UI Manager

```
1 def setup_layout(self):
2     main_layout = QVBoxLayout(self.mw.centralWidget())
3
4     # Set up navigation
5     nav_layout = QHBoxLayout()
6     self.mw.prev_image_btn = QPushButton("Previous Image")
7     self.mw.prev_image_btn.clicked.connect(self.mw.prev_image)
8     self.mw.next_image_btn = QPushButton("Next Image")
9     self.mw.next_image_btn.clicked.connect(self.mw.next_image)
10
11    # Set up content area with image display and control panel
12    content_layout = QHBoxLayout()
13    image_scroll = QScrollArea()
14    image_scroll.setWidgetResizable(True)
15    # [Layout setup continues...]
```

Image Display Manager


```
1 def redraw_image_with_overlays(self):
2     if self.mw.current_image_index < 0 or not self.mw.images:
3         return
4
5     image_data = self.mw.images[self.mw.current_image_index]
6     image = image_data['image'].copy()
7
8     # Convert to OpenCV format for drawing
9     display_cv_img = image.copy()
10
11    # Draw calibration points, center, rings, etc.
12    mc = self.mw.measurement_controller
```

Data Storage Module


```
1 def export_to_csv(self):
2     if not self.measurements:
3         QMessageBox.warning(self, 'Warning', 'No measurements to export')
4         return
5
6     file_path, _ = QFileDialog.getSaveFileName(
7         self,
8         'Save Measurements',
9         '',
10        'CSV Files (*.csv)')
11
12     if not file_path:
13         return
14
15     L = 150
16     n = 1.46
17
18     with open(file_path, 'w', newline='') as f:
19         writer = csv.writer(f, delimiter='\t')
20
21         writer.writerow([
22             'I(A)', 'B(G)',
23             'R_i(mm)', 'R_c(mm)', 'R_o(mm)',
24             'd_i(deg)', 'd_c(deg)', 'd_o(deg)',
25             'B_i(deg)', 'B_c(deg)', 'B_o(deg)',
26             'Δλ_i(nm)', 'Δλ_o(nm)',
27             'ΔE_i(eV)', 'ΔE_o(eV)'
28         ])
29
30     for m in self.measurements:
31         B_field = m.B_field * 1e4
32
33         try:
34             slope, intercept = self.calibration_window.calibration_params
35             current = (B_field - intercept) / slope
36         except (AttributeError, TypeError):
37             QMessageBox.warning(self, 'Warning', 'Please calibrate the magnetic field first')
38         return
39
40
41 def calc_angles(r_mm):
42     if r_mm is None:
43         return None, None
44     alpha = np.arctan(r_mm / L)
45     beta = np.arcsin(np.sin(alpha) / n)
46     return alpha, beta
47
48 alpha_i, beta_i = calc_angles(m.R_inner)
49 alpha_c, beta_c = calc_angles(m.R_center)
50 alpha_o, beta_o = calc_angles(m.R_outer)
51
52 def format_val(val, precision=6):
53     return f'{val:.{precision}f}' if val is not None else ''
54
55 row = [
56     format_val(current, 2),
57     format_val(B_field, 2),
58     format_val(m.R_inner, 3) if m.R_inner else '',
59     format_val(m.R_center, 3) if m.R_center else '',
60     format_val(m.R_outer, 3) if m.R_outer else '',
61     format_val(np.degrees(alpha_i), 4) if alpha_i is not None else '',
62     format_val(np.degrees(alpha_c), 4) if alpha_c is not None else '',
63     format_val(np.degrees(alpha_o), 4) if alpha_o is not None else '',
64     format_val(np.degrees(beta_i), 4) if beta_i is not None else '',
65     format_val(np.degrees(beta_c), 4) if beta_c is not None else '',
66     format_val(np.degrees(beta_o), 4) if beta_o is not None else '',
67     format_val(m.delta_lambda_i * 1e9, 3) if m.delta_lambda_i else '', # nm
68     format_val(m.delta_lambda_o * 1e9, 3) if m.delta_lambda_o else '', # nm
69     format_val(m.delta_E_i / 1.602176634e-19, 6) if m.delta_E_i else '', # eV
70     format_val(m.delta_E_o / 1.602176634e-19, 6) if m.delta_E_o else '' # eV
71 ]
72
73 writer.writerow(row)
74
75 QMessageBox.information(self, 'Success', f'Measurements exported to {file_path}')
```

Usage



In this section, we compare the program's results with theory, verifying the Zeeman effect



Results




Results (continued)





Results (continued)

	R _i (mm)	R _c (mm)	R _o (mm)	Δλ _i (nm)	Δλ _o (nm)	ΔE _i (eV)	ΔE _o (eV)
1	2.042	2.470	2.937	-0.013	0.017	-3.875e-05	5.065e-05
2	1.887	2.451	3.054	-0.016	0.022	-4.910e-05	6.660e-05



Results (continued)



Results (continued)

I(A)	B(G)	R_i(mm)	R_c(mm)	R_o(mm)	α_i(deg)	α_c(deg)	α_o(deg)
β_i(deg)	β_c(deg)	β_o(deg)	Δλ_i(nm)	Δλ_o(nm)	ΔE_i(eV)	ΔE_o(eV)	
0.50	5000.00	5.000	6.000	7.000	1.9092	2.2906	2.6719
1.8297	-0.074	0.087	-0.000220	0.000260	1.3075	1.5687	
1.00	10000.00	5.500	6.500	7.500	2.0999	2.4813	2.8624
1.9601	-0.080	0.094	-0.000240	0.000280	1.4381	1.6992	
1.50	15000.00	6.000	7.000	8.000	2.2906	2.6719	3.0529
2.0905	-0.087	0.100	-0.000260	0.000300	1.5687	1.8297	
2.00	20000.00	6.500	7.500	8.500	2.4813	2.8624	3.2433
2.2208	-0.094	0.107	-0.000280	0.000320	1.6992	1.9601	

Thank you