

Dead time of the GM tube

Slides prepared by:
Abdallah W. Mahmoud


PHYS 3052 – Nuclear Physics Lab

Ahmed Hamed , Hussein
El-Samman





Types of Radiation Detectors


- Method based on the detection of free charge carriers by the ionization, i.e., ionization chamber, proportional counter, Geiger-Muller counter and semiconductor detectors.
- Method based visualization of the tracks of radiation, i.e., Wilson Cloud Chamber, Bubble Chamber, nuclear emulsion plates, spark chamber etc.
- Method based on light sensing, i.e., scintillation counter, cerenkov counter etc





Appetizer






- Sequential edge triggered circuits













GM tube



blue LED





HV supply Converts the 9 V battery voltage to the 400 V needed by the GM tube.

GM tube Detects ionizing radiation: emits a current pulse whenever an ionization event occurs inside the tube.

inverter Converts the current from the GM tube into an inverted voltage pulse.

pulse stretcher Converts the very short pulse from the inverter into a 1.5 ms pulse.

blue LED Flashes every time a pulse occurs.

piezo speaker Produces a click every time a pulse occurs.

filter Accumulates charge from the pulses to produce a voltage roughly proportional to the count rate.




red/green LED Displays a visual indication of the count rate that continuously fades from green to orange to red as the rate increases.

red pump Drives the red LED with a current proportional to the voltage from the filter.

difference amplifier Subtracts the filter voltage from a reference voltage to produce a signal that goes down as the count rate goes up.

green pump Drives the green LED with a current proportional to the voltage from the difference amplifier.






Inside a GM tube

- The Geiger–Müller tube is filled with an inert gas such as helium, neon, or argon at low pressure, to which a high voltage is applied.
- The tube briefly conducts electrical charge when a particle or photon of incident radiation makes the gas conductive by ionization .
- The ionization is considerably amplified within the tube and fed to the processing and display electronics
- When a charged particle or gamma-radiation enters the tube, the argon gas becomes ionized. This triggers a whole avalanche of ions between the electrodes.
- For a brief moment, the gas conducts and a pulse of current flows in the circuit

Pulse stretcher

- The output from the GM tube is very brief – if it were used to drive other circuit components directly, the effect would be nearly imperceptible. The pulse stretcher uses a very useful component known as a 555 timer to turn a short input pulse (of arbitrary shape and duration) into a square pulse.



Dead time in pulse stretcher

- During the time when the pulse is high, the circuit is designed to ignore any subsequent discharges
- This period can be considered a blackout period during which the pulses are processed
- Some pulses may need to be shaped or amplified, so the blackout time can be even bigger
- Some pulse stretchers use a threshold to ignore pulses below a certain voltage




Basic definition

- An instrument that responds sequentially to individual events requires a certain minimum time to recover from one event before it is ready to respond to the next
- This interval, called the dead (or resolving) time, can be due to physical processes in the detector and to instrument electronics
- When a radioactive sample is counted, there is a possibility that two interactions in the detector will occur too close together in time to be registered as separate events

How to model dead time behaviour?

- Two models were proposed
- Following a count, a *paralyzable* detector is unable to provide a second response until a certain dead time τ has passed without another event occurring. Another event during τ causes the insensitive period to be restarted
- A *non paralyzable* detector, on the other hand, simply ignores other events if they occur during τ .



- Such instruments thus actually count the number of intervals between events to which they respond, rather than the number of events themselves
- In practice, counting systems often exhibit behavior intermediate to the two extremes illustrated in the previous slide
- Dead time is one of the main limitations of GM tubes, especially if the 'avalanche' is too overwhelming for them.

Calculations

- For r_t true events, the probability of the GM tube being ready to record is $1 - r_t \tau$
- This probability is equal to the fraction of the detected count rate r_c / r_t
- We equate both sides, getting $r_c / r_t = 1 - r_t \tau$
- Rearranging the equation, we get the expression on the right for a non paralyzable detector

$$r_t = \frac{r_c}{1 - r_c \tau} \quad (\text{Nonparalyzable}).$$

$$(r_c \tau \ll 1),$$

$$r_t \cong r_c (1 + r_c \tau).$$

A paralyzable detector is different...

- Only intervals longer than τ are registered
- We need the distribution of time intervals between successive random events that occur at the average rate r_t
- The average number of events that take place in a time t is $r_t \cdot t$
- If an event occurs at time $t = 0$, then the probability that no events occur in time t immediately following that event is given by the Poisson term, $P_0 = \exp(-r_t \cdot t)$
- The probability that an event will occur in the next time interval dt is $r_t dt$
- Therefore, given an event at time $t = 0$, the probability that the next event will occur between t and $t + dt$ is:

$$P(t) dt = r_t e^{-r_t t} dt.$$

Calculations


- The probability that a time interval larger than t will elapse is $\rightarrow\rightarrow\rightarrow$
- This probability is equal to the fraction of the detected count rate r_c/r_t
- We equate both sides, getting $r_c/r_t = \exp(-r_t t)$
- Rearranging the equation, we get the expression on the right for a paralyzable detector

$$\int_{\tau}^{\infty} r_t e^{-r_t t} dt = e^{-r_t \tau}.$$

$$r_c = r_t e^{-r_t \tau} \quad (\text{Paralyzable}).$$

$$(r_t \tau \ll 1),$$

$$r_c = r_t (1 - r_t \tau).$$




More Calculations

- For the nonparalyzable model, as r_t increases, r_c approaches the asymptotic value $1/\tau$
- For the paralyzable model, the observed count rate goes through a maximum $1/e\tau$ when $r_t = 1/\tau$
- With increasing event rates, the measured count rate with a paralyzable system will decrease beyond this maximum and approach zero, because of the decreasing opportunity to recover between events.

$$r_t = \frac{r_c}{1 - r_c \tau} \quad (\text{Nonparalyzable}).$$

$$r_c = r_t e^{-r_t \tau} \quad (\text{Paralyzable}).$$



Example

- A counting system has a dead time of $1.7 \mu\text{s}$. If a count rate of $9 \times 10^4 \text{ s}^{-1}$ is observed, what is the true event rate if the counter is (a) nonparalyzable or (b) paralyzable?
- (a)

$$r_t = \frac{9 \times 10^4 \text{ s}^{-1}}{1 - (9 \times 10^4 \text{ s}^{-1})(1.7 \times 10^{-6} \text{ s})} = 1.06 \times 10^5 \text{ s}^{-1}.$$

- (b)


$$9 \times 10^4 = r_t e^{-1.7 \times 10^{-6} r_t} \quad \Rightarrow \quad 1.7 \times 10^{-6} r_t = \ln \frac{r_t}{9 \times 10^4}.$$

Solving by iteration

$$1.7 \times 10^{-6} r_t = \ln \frac{r_t}{9 \times 10^4}.$$

- We try the solution $r_t = 1.06 \times 10^5 \text{ s}^{-1}$ here. The lefthand side then has the value 0.180, compared with the smaller value on the right-hand side, 0.164
- Since the paralyzable counter misses more events than the nonparalyzable, r_t should be larger now than in part (a)
- Trying $r_t = 1.10 \times 10^5 \text{ s}^{-1}$ gives 0.187 on the LHS compared with the larger 0.201 on the right. Thus, equality is satisfied by a value of r_t between these two trial values. Further refinement yields the solution $r_t = 1.08 \times 10^5 \text{ s}^{-1}$
- Since this solution is close to the result for the nonparalyzable counter, we expect the second solution to be at a higher event rate.
- For $r_t = 10^{-6} \text{ s}^{-1}$, the LHS and RHS give, respectively, 1.70 and 2.41. For $r_t = 10^7 \text{ s}^{-1}$, the results are 17.0 and 4.71. Therefore, the solution is between these two values of r_t . One finds $r_t = 1.75 \times 10^6 \text{ s}^{-1}$.

PULSE AMPLITUDE



- Geiger-Muller tubes exhibit Dead Time effects due to the recombination time of the internal gas ions after the occurrence of an ionizing event
- The counter discharge occurs very close to the wire, and the negative particles, usually electrons, are collected very rapidly
- The positive ions move relatively slowly, so that as the discharge proceeds a positively charged sheath forms around the wire. This has the effect of reducing the field around the wire to a value below that corresponding to the threshold voltage, and the discharge ceases
- If another ionizing event triggers the counter at this stage, a pulse smaller than normal is obtained, as the full voltage across the counter is not operative.
- However, if the positive ions reach the cathode before the next particle arrives, the pulse will be of full size.



Dead time vs Recovery time

- Physically, the deadtime is the time required for the positive ions to travel far enough from the wire to permit the resumption of Geiger action, and the recovery time is the time required for the positive ions to reach the counter wall
- The resolving time is the addition of dead time and recovery time
- During dead time, GM tube is completely unresponsive
- During recovery time, GM tube is partially responsive, and smaller pulses are generated that may or may not be detected
- Resolving time is the minimum time interval by which two pulses must be separated to be detected as separate pulses by the counter

Experimental Procedure

The two-source method





Since $N_s = N_1 + N_2$

$$N_s = \frac{n_1}{(1 - n_1 t)} + \frac{n_2}{(1 - n_2 t)} \quad (5)$$

From (4) we have

$$n_s = \frac{N_s}{(1 + N_s t)} \quad (6)$$

Substituting (5) into (6), we obtain after manipulating:

$$n_s = \frac{n_1 + n_2 - 2n_1 n_2 t}{(1 - n_1 n_2 t^2)} \quad (7)$$

Normally $n_1 n_2 t^2 \ll 1$, we can approximate eq (7) to

$$n_s = n_1 + n_2 - 2n_1 n_2 t \quad (8)$$


which yields

$$t = \frac{n_1 + n_2 - n_s}{2n_1 n_2} \quad (9)$$

$$N_1 = \frac{n_1}{(1 - n_1 t)} \quad (2)$$

$$N_2 = \frac{n_2}{(1 - n_2 t)} \quad (3)$$

$$N_s = \frac{n_s}{(1 - n_s t)} \quad (4)$$



Negative dead time

- Primarily due to errors in the measurement process
- The radioactive source is not properly aligned with the detector, leading to inconsistent count rates
- If count rate not corrected for background, you may end up with a negative value for the dead time
- Weak sources have relatively larger fluctuations in count rates, which can lead to an inaccurate dead time
- Faulty or malfunctioning electronic components could produce inconsistent pulses, leading again to an inconsistent count rate

Some improvements

- Add a halogen gas or organic vapour that increases the speed of the quenching mechanism and stops the ionization avalanche faster
- Use an external quenching circuit to implement electronic quenching
- Due to the dependence of dead time on the operating voltage, choose your operating voltage wisely, preferably one in the middle of the plateau
- Use advanced electronics to optimize the processing of the signal, such that the generated pulse can be narrower, reducing the electronic dead time
- Optimize the geometry of the GM tube such the anode is not completely covered by the positively charged cloud
- Make sure the source has a robust setup such that the alignment between the source and detector is unchanged throughout the experiment

Thank you

